



QUICK-START GUIDE




Start-up Guide
IMPS-IMVS Optical Bench

1

Install software as described in other Gamry *Quick-Start Guides*, then shut down computer.

2

Assemble two posts for cage supports.



Adjust height of cage support, then tighten thumbscrew on post holder.

3

Assemble post for photodiode.

Photodiode

Lens tube

Lens-tube
slip ring

Optical post

Post holder

Rail carrier

1/4"-20
cap screw



Adjust position of
photodiode, then
tighten locking
screw.

Adjust height of
cage support,
then tighten
thumbscrew on
post holder.



4

Assemble post for LED.

Lens-tube
slip ring

LED assembly

Optical post

Post holder

Rail carrier

1/4"-20
cap screw




Adjust position of
LED, then tighten
locking screw.

Adjust height of
cage support,
then tighten
thumbscrew on
post holder.




5 Assemble base plate.




6

Set up posts on base plate in this order.



Equalize heights of LED and photodiode using scales on optical posts. Align LED and photodiode to center of cage posts.


Affix all posts by tightening thumbscrews on rail carriers.




7 Assemble cage system.

Connect 2 × short rods together (2 sets).
Slide 2 × rod assemblies through upper holes.

Slide 2 × long rods
through lower holes.




Affix all rods by
tightening 8 ×
locking screws
on each cage
plate.




8

Attach BNC cable to photodiode.



9

Attach breakout cables to ends of both LED and photodiode cables.



Follow color-conventions as shown.

10

Connect photodiode breakout cables to Serf potentiostat.



Follow color-conventions as shown.

11

Connect LED breakout cables to Master potentiostat.



Make sure the potentiostats are OFF. Turn on computer. Switch the Main Power of the IPH ON. Microsoft Windows® detects your IPH and dongle, and a **Found New Hardware** prompt appears. If necessary, select **Install Software Automatically**.

12 When setup is complete, double-click the **Gamry Framework** icon. A **Device Settings** window appears for the USB dongle.




13 On the packing slip, match the dongle's serial number to the proper section for software authorization codes. Enter these codes by clicking **Add** and typing the package name and 10-digit authorization code for each package.



14 Turn on your first potentiostat. Your computer recognizes the instrument. Click **Install Software Automatically** if such a dialog box appears. When presented with the **Device Settings** dialog box, change the potentiostat name if desired. Do not install authorization codes for the potentiostats that are part of the MultiEchem.

After a moment, your potentiostat should appear next to the **Devices Present** along with a green virtual LED. Repeat for additional potentiostats.





15 Follow the steps in *Quick-Start Guide #2: USB Potentiostat Calibration* to calibrate your potentiostats. Calibrate the LED.

We have a variety of resources available to help you get started. Feel free to visit our website to find out more information on:

Application Notes - <http://www.gamry.com/application-notes/>

Technical Support - <http://www.gamry.com/service-support/>

Training Videos - <http://www.youtube.com/gamryinstruments>

WHAT DOES GAMRY SOFTWARE DO?



Gamry Framework™

Controls your Gamry potentiostat for advanced and flexible data acquisition. Select from standardized experiments that are grouped by research type or use the Sequence Wizard to build complex automated experiments.



Echem Analyst™

Quick and easy data analysis. Open your data files with Echem Analyst to perform specialized analysis algorithms and produce high quality plots. Plots can be customized, overlaid, and scaled, or data can be exported for use in other plotting programs.



My Gamry Data™

The default data folder location for Gamry Framework data files. You will find a shortcut on your desktop after installation. The data folder location can be changed within Gamry Framework by selecting: Options > Path.