



QUICK-START GUIDE



Guide #1 Multichannel System Installation

1 Insert installation media and click Install Software.




Be sure to have Administrative rights on your computer before installation. If you do not, contact your IT department for help.

2 The Gamry Software Installation program runs.

NOTE: If you have Gamry software previously installed, you are asked to remove previous versions of the software and the Gamry device drivers. Click **YES**; all previous data will be saved.

- When asked to select folder location, click **Next**.
- Follow prompts to finish installation. Restart computer.

3 Open the Gamry Framework™. Gamry Instrument Manager software automatically opens, showing the new instrument and its characteristics.



4

To change the label of an instrument, click the pencil icon next to the instrument's label.




5

Close the **Instrument Manager**.

6

After a moment, your potentiostat should appear next to **Devices Present** along with a green virtual LED. Repeat for additional potentiostats.



7

Next, follow the steps in the *Quick-start Guide #2: USB Potentiostat Calibration* to calibrate your potentiostats.

Be sure to check our website,
www.gamry.com/support/software-updates/
for the most current updates to your software.



WHAT DOES GAMRY SOFTWARE DO?



Gamry Framework™

Potentiostat control for flexible data acquisition. Select from standardized experiments grouped by research type, or use the Sequence Wizard to build complex automated experiments.



Echem Analyst™

Quick and easy data analysis. Open data files with Echem Analyst for specialized analysis algorithms and high-quality plots. Customize, overlay, and scale plots, or export data.



My Gamry Data™

The default data-folder location for Gamry Framework, with a shortcut on your desktop after installation. Change the folder location within Gamry Framework via **Options > Path**.



Virtual Front Panel™

Software-based front panel for quick access to Gamry potentiostats' functions, like a front panel of an early analog potentiostat; and to perform simple electrochemical experiments.



Electrochemical Signal Analyzer™

Designed specifically for the acquisition and analysis of time-dependent electrochemical noise signals.



Resonator™

Data-acquisition and -control software for the Gamry eQCM™. Contains a full suite of physical electrochemistry techniques.



Electrochemistry Toolkit™

A sophisticated package for complete access to the capabilities of Gamry potentiostats in the software environment of your choice.