


Electrochemical Impedance Spectroscopy: Quick Check of System Accuracy

Gamry Instruments



EIS: Accuracy Contour Plot vs. Quick Check

- A complete ACP shows total, normal system performance (see right)
 - Generally collected using normal operating parameters (e.g. 10 mV potentiostatic EIS)
 - May also be optimized for each point/region
 - Takes a long time to measure and determine
 - (Open and Shorted Lead EIS spectra are *Not* the ACP)
- EIS Quick Check:
 - Simple test that will work across platforms
 - Gives a rough check of system performance
 - Establishes proper setup for making accurate measurements



How to Run an EIS Quick Check

- Run EIS on Some Precision Resistors
 - Potentiostatic, 10 mV AC
 - 100 Ω 1 k Ω and 10 k Ω
 - Shows Real System Response
- Run Shorted Lead EIS
 - Galvanostatic, 100-400 mA AC
 - Shows Low Impedance Ability
 - Measure of Cable Inductance
- [Optional] Run an Open Lead EIS
 - Potentiostatic, 50 mV AC
 - Measure of Cable Capacitance, System Limit
 - Not Necessarily Useful for Accuracy


Cable Setup Matters

- For some impedance measures Cable positioning makes a big difference.
- 1 k Ω : Any should work
- 100 Ω and Shorted: minimize inductance by twisting cables (current pair together from one side, sense pair together from the other)
- 10 k Ω and Open: minimize capacitance by separating work leads from counter/ref leads as much as possible
 - (for Open connect work to work sense and counter to reference, place working pair inside Faraday cage)




Good Resistor Response

- 100 Ω : Green
- 1 k Ω : Blue
- 10 k Ω : Red
- Magnitude w/in 1%
- Phase w/in 2°
 - Down capacitive
 - Up Inductive




Bad Resistor Response

- Gamry: Blue
- *Other*: Red
- *Other* 1 MHz System
 - 40% Error in Z
 - 30° Phase Error
 - System only accurate to ~50 kHz




Shorted Lead Curve

- Very small resistor
- Inductance will be present at higher f
- Note:
 - Calculated resistance of copper braid is $10 \mu\Omega$ per mm, with 2-3 mm between leads (or about the measured $30 \mu\Omega$)



Open Lead Curve

- Should look like a capacitor
 - 90° phase
 - Capacitance based on cabling
 - Low f rolloff may be present on less sensitive systems



Quick Check Take Home

- Instruments have limits
- Setup matters
- Know where your system lies
 - High Impedance (Coatings)
 - Low Impedance (Energy Devices)
 - Other...(Corrosion, Sensors)